

## Annexure - VI

### List of instruments with make and model

Sr. No	Name of instrument	Make	Model	INST.ID. NO	Location (Floor/ Room No.)
1	HPLC	Shimadzu	SPD-M10ATVP	CDTL/LC/01	F-117
2	HPLC	Shimadzu	Prominance	CDTL/LC/02	F-117
3	HPLC	Merck Hitachi	L7400	CDTL/LC/03	F-117
4	HPLC	Merck Hitachi	L7400	CDTL/LC/04	F-117
5	HPLC	Agilent	1100	CDTL/LC/05	F-125
6	HPLC	Waters	2996	CDTL/LC/06	F-125
7	HPLC	Agilent	1200	CDTL/LC/07	F-125
8	HPLC	Agilent	1200	CDTL/LC/08	F-125
9	HPLC	Perkin Elmer	Flexar	CDTL/LC/09	G-013
10	HPLC	Perkin Elmer	Flexar	CDTL/LC/10	G-013
11	HPLC	Perkin Elmer	Flexar	CDTL/LC/11	G-013
12	HPLC	Perkin Elmer	Flexar	CDTL/LC/12	G-013
13	HPLC	Dionex	3000	CDTL/LC/13	G-013
14	HPLC	Dionex	Ultimate-3000	CDTL/LC/14	G-013
15	HPLC	Dionex	Ultimate-3000	CDTL/LC/15	G-013
16	HPLC	Dionex	Ultimate-3000	CDTL/LC/18	S-203
17	HPLC	Dionex	Ultimate-3000	CDTL/LC/19	S-203
18	HPLC	Dionex	Ultimate-3000	CDTL/LC/20	S-203
19	HPLC	Dionex	Ultimate-3000	CDTL/LC/21	T-302
20	HPLC	Dionex	Ultimate-3000	CDTL/LC/22	T-302
21	HPLC	Dionex	Ultimate-3000	CDTL/LC/23	T-302
22	HPLC	Agilent	1260	CDTL/LC/17	G-013
23	Fast HPLC (UPLC)	Dionex	Ultimate-3000	CDTL/LC/16	G-013
24	GC	ThermoFinnigan	Trace	CDTL/Chem/GC/01	F-121
25	GC	Agilent	7890A	CDTL/Chem/GC/02	F-121
26	GC	ThermoFinnigan	Trace	CDTL/Chem/GC/03	G-013
27	HPTLC	Camag	Linomat-IV	CDTL/HPTLC/01	F-120
28	HPTLC	Camag	Linomat-V	CDTL/HPTLC/02	F-120
29	TLC(plate Heater)	Camag	NA	CDTL/TLC/01	F-120
30	TLC-Dryer	Camag	NA	CDTL/DR/01	F-121
31	Scanner	Camag	NA	CDTL/HPTLC-SCN/01	F-122
32	Visulizer	Camag	NA	CDTL/HPTLC-VIS/01	F-123
33	TLC-UV Cabinet	Camag	NA	CDTL/UVC/01	F-124
34	TLC-UV Cabinet	Camag	NA	CDTL/UVC/01	F-125
35	Autotitrator	Metrohm	877 Titriplus(702)	CDTL/AUTO/01	F-104
36	Autotitrator	Metrohm	888 Titrand	CDTL/AUTO/03	GF-01
37	Autotitrator	Metrohm	888 Titrand	CDTL/AUTO/04	F-102
38	Autotitrator	Labindia	Titra smart	CDTL/AUTO/05	GF-01
39	Autotitrator	Labindia	Titra smart	CDTL/AUTO/06	GF-01
40	KF	Metrohm	795KFT	CDTL/KF/01	F-106
41	KF	Metrohm	890KFT	CDTL/KF/02	F-106
42	KF	Labndia	KAFI smart	CDTL/KF/03	GF-01
43	Coulometer	Metrohm	831KFC	CDTL/Coul/01	F-106
44	Atomic Absorption Spectrophotometer	GBC	Savant AA	CDTL/AAS/2	F-122
45	Autoclave	Osworld	NA	CDTL/MIC/AUTO/01	S-205
46	Autoclave	Medica Lab Equipment	325.91.6	CDTL/MIC/AUTO/02	S-225
47	Autoclave	Osworld	NA	CDTL/MIC/AUTO/03	S-225
48	Autoclave	Preccision	PB-08	CDTL/AUTOCLAVE/04	S-225
49	Analytical Balance	LCGC Radwag	AS220/C/2	CDTL/BAL/04	F-104
50	Analytical Balance	LCGC Radwag	AS220/C/2	CDTL/BAL/06	F102
51	Analytical Balance	LCGC	377317/12	CDTL/BAL/07	S-202
52	Analytical Balance	Mettler toledo	AB204	CDTL/BAL/08	S-203
53	Analytical Balance	Mettler toledo	AT200	CDTL/BAL/09	S-204
54	Analytical Balance	Sartorius	BSA224S-CW	CDTL/BAL/11	S-224
55	Analytical Balance	Denver	SI-203	CDTL/BAL/10	S-220
56	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/15	F-128
57	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/16	F-106
58	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/17	S-204
59	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/18	G-02

60	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/19	S-204
61	Balance	Citizen	-----	CDTL/BAL/20	S-203
62	Balance	Citizen	-----	CDTL/BAL/21	F-102
63	Balance	Sartorius	Secura12510IN	CDTL/BAL/22	F-108
64	Balance	Sartorius	Secura12510IN	CDTL/BAL/23	F-108
65	BET Kinetic instrument	Lonza	ELX-808	CDTL/KIN/01	S-217
66	BOD Incubator	Osworld	3979	CDTL/BOD/01	S-223
67	BOD Incubator	Osworld	3980	CDTL/BOD/02	S-223
68	BOD Incubator	Metalab	JRIC-10	CDTL/BOD/03	S- 223
69	BOD Incubator	REMI	CI-6S LED GMP	CDTL/BOD/04	S- 223
70	Biosafety Cabinet	-----	-----	CDTL/BSC/01	S-221
71	Break load test apparatus	Lab fine	-----	CDTL/COS-BLT/01	T-303
72	Centrifuge	Electrolab	CM-6MT	CDTL/CENTRI/01	F-101
73	Centrifuge	Remi motors	R-24	CDTL/CENTRI/02	G-010
74	Centrifuge	Remi motors	R-8C BL	CDTL/CENTRI/03	G-03
75	Centrifuge	Remi motors	R-8M PLUS	CDTL/CENTRI/04	F-105
76	Centrifuge	Remi motors	R-8C PLUS	CDTL/CENTRI/05	S-202
77	Centrifuge	Remi motors	R-8C PLUS	CDTL/CENTRI/06	T-304
78	Compressor	NA	NA	CDTL/Compsr/01	T-303
79	Dehumidifier	Sharp	NA	CDTL/DH/03	G-001
80	Dehumidifier	Sharp	NA	CDTL/DH/04	F-104
81	Disintegration apparatus	Veego	VTD-AV	CDTL/DT/02	S-202
82	Disintegration apparatus	Electrolab	EDI-2SA	CDTL/DT/05	F-119
83	Disintegration apparatus	Electrolab	EDI-2SA	CDTL/DT/06	T-304
84	Digital Force Gauge	Chattillon	DFA-010	CDTL/DFG/01	T-303
85	Digital Screw Gauge	Mitutoyo	293-561-30	CDTL/DSG/01	T-303
86	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/01	F-119
87	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/02	F-119
88	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/03	S-205
89	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/04	G-04
90	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/05	G-04
91	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/06	G-04
92	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/07	T-305
93	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/08	G-04
94	Dissolution apparatus	LabIndia	DS8001	CDTL/Disso/09	S-204
95	Dissolution media Prep sys	Lab inida	Media MET+	CDTL/DMP/01	G-04
96	Digital Flow meter for HPLC	GJC Instrument	5025000	CDTL/FM/01	F-125
97	Distillation unit(Food lab)	Velp Scientific	NA	NA	G-10
98	Drying time Machine	Labfine	NA	NA	T-304
99	Flame photometer	NA	-----	CDTL/FP/01	F-122
100	Foam Height Apparatus	NA	NA	CDTL/FHA/01	S-213
101	Freezer (Ultra deep)	SR lab Instruments	-----	CDTL/RF/039	S-217
102	FTIR	Shimadzu	IR Affinity - 1	CDTL/FTIR/02	F-118
103	FTIR-Press	NA	-----	CDTL/RP/01	118
104	FTIR-Press	NA	-----	CDTL//RP/02	118
105	FTIR-Press	Technosearch	-----	CDTL//RP/03	118
106	Fume Hood	Labguard	ST1500ABP	CDTL/FH/01	F-106
107	Fume Hood	Labguard	ST1500ABP	CDTL/FH/02	S-205
108	Fume Hood	Labguard	ST1500ABP	CDTL/FH/03	S-205
109	Fume Hood	Prama tech	NA	CDTL/FH/04	T-305
110	Funnel shaker for solvent extractor	Jiotech	ST1500ABP	CDTL/SFS/01	F-117
111	Gel Electrophoresis	Himedia	NA	CDTL/GEL/01	S-202
112	Heating mantle	NA	NA	CDTL/HM/03	S-225
113	Heating mantle	NA	NA	CDTL/HM/04	F-107
114	Heating mantle	jindal	NA	CDTL/HM/05	S-205
115	Heating mantle	NA	NA	CDTL/HM/06	G-10
116	Heating mantle	Labquest	HME500	CDTL/HM/07	F-117
117	Heating mantle	Labquest	HME500	CDTL/HM/08	F-117
118	Heating mantle	Labquest	HME500	CDTL/HM/09	S-205
119	Heating mantle	Labquest	HME500	CDTL/HM/10	T-305
120	Heating Block	-----	NA	NA	S-217
121	Heating Block	-----	NA	NA	S-217
122	Hot plate digital	Tarson	NA	CDTL/HP/03	F-107
123	Hot plate digital	Tarson	NA	CDTL/HP/04	S-207

124	Incubator	Pathak		CDTL/INC/01	S-223
125	Incubator	Metlab	MSI-67	CDTL/INC-02	S-223
126	Incubator	Metlab	MSI-67	CDTL/INC-03	S-223
127	Laminar Air Flow	Air borne	NA	CDTL/LAF/01	S-219
128	Laminar Air Flow	Air borne	NA	CDTL/LAF/02	S-220
129	Laminar Air Flow	Air borne	NA	CDTL/LAF/03	S-221
130	Laminar Air Flow	Air borne	NA	CDTL/LAF/04	S-222
131	Magnetic stirrer	Emtec	NA	CDTL/MAG/02	S-202
132	Magnetic stirrer	Remi	NA	CDTL/MAG/03	G-04
133	Magnetic stirrer	Labquest	NA	CDTL/MAG/04	F-104
134	Magnetic stirrer	Labquest	NA	CDTL/MAG/05	T-305
135	Magnetic stirrer	Labquest	NA	CDTL/MAG/06	S-205
136	Melting point apparatus	LabIndia	MR-VIS	CDTL/MP/01	F-124
137	Microscope	Labomed	NA	CDTL/MIC/01	S-217
138	Microscope	Labomed	Vision2000	CDTL/MIC/02	S-217
139	Microscope	Labomed	Vision2001	CDTL/MIC/03	S-218
140	Mixer	Kay	NA	CDTL/MIX/02	T-305
141	Muffle furnace	NA	NA	CDTL/MF/01	205
142	Muffle furnace	NA	NA	CDTL/MF/02	205
143	Needle burner & Syringe destroyer	Hi Media	NA	CDTL/NSD/01	
144	Oil Bath	Metalab	NA	CDTL/OB/01	T-304
145	Oven	PATHAK	NA	CDTL/OV/02	F-107
146	Oven	METALAB	NA	CDTL/OV/01	F-101
147	Oven	JINDAL	NA	CDTL/OV/09	G-03
148	Oven	JINDAL	NA	CDTL/OV/10	G-04
149	Oven	MIC	NA	CDTL/OV/04	F-102
150	Oven	MIC	NA	CDTL/OV/08	S-202
151	Oven	Pathak	NA	CDTL/OV/05	S-205
152	Oven	Pathak	NA	CDTL/OV/06	S-224
153	Oven	MIC	NA	CDTL/OV/07	S-224
154	Oven	i therm	NA	CDTL/OV/03	S-205
155	Orbital shaker	Remi	NA	CDTL/SH/01	S-203
156	Orbital shaker	Remi	NA	CDTL/SH/02	G-10
157	Payoff measuring Tester	Labfine	NA	CDTL/PO/01	
158	pH Meter	LabIndia	Pico Plus	CDTL/pH/03	F-119
159	pH Meter	Euetech inst	pH 510	CDTL/pH/04	S-224
160	pH Meter	Euetech inst	Ph 700	CDTL/pH/09	F-104
161	pH Meter	ANALAB	pHcal10	CDTL/pH/12	F-104
162	pH Meter	Agilent	3200P	CDTL/pH/10	G-04
163	pH Meter	Agilent	3200P	CDTL/pH/11	S-202
164	pH Meter	Horiba	F72A	CDTL/pH/13	F-119
165	pH Meter	-----	-----	CDTL/pH/14	S- 201
166	Polarimeter	Rudolph research	-----	CDTL/POL/01	F-117
167	Polarimeter	Anton-Par	-----	CDTL/POL/02	F-117
168	Polarimeter	Anton-Par	-----	CDTL/POL/03	S-205
169	Profile projector	SIPCON	SP-300	CDTL/PP/02	T-304
170	Particle counter	-----	-----	CDTL/PC/01	S-203
171	Peel off Tester	Lab fine	-----	CDTL/COS-PT/01	T-303
172	Refractometer	Atago	RX-7000i	CDTL/RI/02	F-117
173	Scratch hardness tester	Labfine	NA	CDTL/SHT/01	ASN/Ramesh
174	Softening point apparatus	Labfine	NA	CDTL/SPA/01	ASN/Ramesh
175	Sonicattor	Pci service	NA	CDTL/Soni/02	F-104
176	Sonicattor	MIC	-----	CDTL//Soni/04	T-304
177	Sonicattor	MIC	-----	CDTL//Soni/05	S-203
178	Sonicattor	Labdevice	NA	CDTL//Soni/07	F101
179	Sonicattor	MIC	-----	CDTL//Soni/08	S-202
180	Sonicattor	MIC	-----	CDTL//Soni/09	S-204
181	Sonicattor	Pci service	NA	CDTL//Soni/10	G-03
182	Sonicattor	Pci service	NA	CDTL//Soni/11	F-101
183	Solvent cabinet	Prama	-----	-----	G-011
184	Solvent cabinet	Prama	-----	-----	G-011
185	Spot Extractor	Prama	-----	CDTL/SPEX/01	F-104
186	Spot Extractor	Prama	-----	CDTL/SPEX/02	F-104
187	Spot Extractor	Prama	-----	CDTL/SPEX/03	F-101

188	Spot Extractor	Prama	-----	CDTL/SPEX/04	F-101
189	Spot Extractor	Prama	-----	CDTL/SPEX/05	S-202
190	Spot Extractor	Prama	-----	CDTL/SPEX/06	S-203
191	Spot Extractor	Prama	-----	CDTL/SPEX/07	S-204
192	Standard Weight Box	IDEMI Mumbai	NA	-----	QA
193	Tensile compression Device	Mecmesin	AFG-100N	CDTL/TCD/03	T-303
194	Traceble Digital Calliper	Fisher Scientific	NA	CDTL/TDC/01	T-303
195	Traceble Digital Timer	Fisher Scientific	NA	CDTL/TDT/01	T-303
196	TOC Analyzer	shimazu	TOC-LCPH	CDTL/TOC/1	G-02
197	Tablet Friability Tester	lab india	-----	CDTL/TFT/01	G-013
198	Tablet Hardnes Tester	lab india	-----	CDTL/THT/01	G-013
199	UV/VIS SPECTROPHOTOMETER	LabIndia	UV300+	CDTL/UV/01	G-03
200	UV/VIS SPECTROPHOTOMETER	LabIndia	UV300+	CDTL/UV/02	S-203
201	UV/VIS SPECTROPHOTOMETER	Perkin Elmer	UV150-02	CDTL/UV/03	F-117
202	DOUBLE BEAM UV/VIS SPECTROPHOTOMETER	Agilent	Cary 60	CDTL/UV/04	T-302
203	Ultra violet Nail paint Dryer	-----	-----	CDTL/COS-ND/01	T-303
204	Vaccume oven	Osworld	NA	CDTL/VO/01	F-105
205	Vaccume oven	LAB HOSP	NA	CDTL/VO/02	F-107
206	Vaccume oven	HMG India	NA	CDTL/VO/03	F-105
207	Vaccume oven	HMG	NA	CDTL/VO/04	S-203
208	Vaccume oven	Fourtech	NA	CDTL/VO/05	S-206
209	Vaccume oven	SM Scientific	NA	CDTL/VO/06	F-105
210	Vaccume oven	META LAB	NA	CDTL/VO/07	G-02
211	Vaccume oven	JEIOTECH	OV-12	CDTL/VO/08	F-105
212	Vortex	SPINIX	NA	CDTL/VOR/01	F-102
213	Vortex	SPINIX	NA	CDTL/VOR/02	F-202
214	Vortex	Remi	NA	CDTL/VOR/03	S-204
215	Vortex	SPINIX	NA	CDTL/VOR/04	S-220
216	Vortex	Remi	NA	CDTL/VOR/05	S-220
217	Vortex	-----	-----	CDTL/VOR/07	S-217
218	Vortex	Remi	CM101	CDTL/VOR/08	F104
219	Vortex	Remi	CM102	CDTL/VOR/09	G-03
220	Vortex	Labquest	VMSP	CDTL/VOR/10	F104
221	Vortex	Labquest	VMSP	CDTL/VOR/11	T-305
222	Vortex	Labquest	VMSP	CDTL/VOR/12	S-205
223	Vortex	Labquest	VMSP	CDTL/VOR/13	G-10
224	Vortex mix	Remi	NA	CDTL/VOR/14	G-04
225	Water Bath-Shaker	-----	NA	CDTL/SWB/01	S-202
226	Water Bath-Shaker	-----	NA	CDTL/SWB/02	F-107
227	Water Bath	-----	NA	CDTL/WB/06	S-224
228	Water Bath	-----	NA	CDTL/WB/10	S-205
229	Water Bath	-----	NA	CDTL/WB/11	F-107
230	Water Bath	Foodlab	NA	CDTL/WB/12	F-107
231	Water Bath	Foodlab	NA	CDTL/WB/13	T-305
232	Water Bath	Labquest	WBC-012	CDTL/WB/14	F-107
233	Water Bath	Labquest	WBC-012	CDTL/WB/15	S-205
234	Water Bath	Labquest	WBC-012	CDTL/WB/16	T-305
235	Water Bath	Labquest	WBC-012	CDTL/WB/17	G-10
236	Water purifier	Pure Lab Ultra	ELGA DV-35	CDTL/WP/01	T-304
237	Water purifier	Millipore	Q815 V	CDTL/WP/02	F-101
238	Water purifier	Millipore	Elix	CDTL/WP/05	F-105
239	Water purifier	Thermo	Genpure	CDTL/WP/04	G-13
240	Water purifier	Thermo	Genpure	CDTL/WP/03	S-211
241	Zone Reader	-----	-----	CDTL//ZR/01	S-227
242	Rotational Viscometer	Anton-Paar	ViscoQC 300	CDTL/VM/02	T-303
243	Density Meter	Anton-Paar	DMA501	CDTL/DM/01	F-117
244	Circulating Water Bath	CENTROFIX'	3S-HCB	CDTL/CWB/01	S-205
245	Chiller Sonicator Bath	CENTROFIX'	SM 200 CH US	CDTL/CSB/01	S-203
246	Chiller Sonicator Bath	CENTROFIX'	SM 200 CH US	CDTL/CSB/02	F-106