

Annexure - VI

List of instruments with make and model

Sr. No	Name of instrument	Make	Model	INST.ID. NO	Location (Floor/ Room No.)
1	HPLC	Shimadzu	SPD-M10ATVP	CDTL/LC/01	F-117
2	HPLC	Shimadzu	Prominance	CDTL/LC/02	F-117
3	HPLC	Merck Hitachi	L7400	CDTL/LC/03	F-117
4	HPLC	Merck Hitachi	L7400	CDTL/LC/04	F-117
5	HPLC	Agilent	1100	CDTL/LC/05	F-125
6	HPLC	Waters	2996	CDTL/LC/06	F-125
7	HPLC	Agilent	1200	CDTL/LC/07	F-125
8	HPLC	Agilent	1200	CDTL/LC/08	F-125
9	HPLC	Perkin Elmer	Flexar	CDTL/LC/09	G-013
10	HPLC	Perkin Elmer	Flexar	CDTL/LC/10	G-013
11	HPLC	Perkin Elmer	Flexar	CDTL/LC/11	G-013
12	HPLC	Perkin Elmer	Flexar	CDTL/LC/12	G-013
13	HPLC	Dionex	3000	CDTL/LC/13	G-013
14	HPLC	Dionex	Ultimate-3000	CDTL/LC/14	G-013
15	HPLC	Dionex	Ultimate-3000	CDTL/LC/15	G-013
16	HPLC	Dionex	Ultimate-3000	CDTL/LC/18	S-203
17	HPLC	Dionex	Ultimate-3000	CDTL/LC/19	S-203
18	HPLC	Dionex	Ultimate-3000	CDTL/LC/20	S-203
19	HPLC	Dionex	Ultimate-3000	CDTL/LC/21	T-302
20	HPLC	Dionex	Ultimate-3000	CDTL/LC/22	T-302
21	HPLC	Dionex	Ultimate-3000	CDTL/LC/23	T-302
22	HPLC	Agilent	1260	CDTL/LC/17	G-013
23	Fast HPLC (UPLC)	Dionex	Ultimate-3000	CDTL/LC/16	G-013
24	GC	Thermofinnigan	Trace	CDTL/Chem/GC/01	F-121
25	GC	Agilent	7890A	CDTL/Chem/GC/02	F-121
26	GC	Thermofinnigan	Trace	CDTL/Chem/GC/03	G-013
27	HPTLC	Camag	Linomat-IV	CDTL/HPTLC/01	F-120
28	HPTLC	Camag	Linomat-V	CDTL/HPTLC/02	F-120
29	TLC(plate Heater)	Camag	NA	CDTL/TLC/01	F-120
30	TLC-Dryer	Camag	NA	CDTL/DR/01	F-121
31	Scanner	Camag	NA	CDTL/HPTLC-SCN/01	F-122
32	Visulizer	Camag	NA	CDTL/HPTLC-VIS/01	F-123
33	TLC-UV Cabinet	Camag	NA	CDTL/UVC/01	F-124
34	TLC-UV Cabinet	Camag	NA	CDTL/UVC/01	F-125
35	Autotitrator	Metrohm	877 Titripplus(702)	CDTL/AUTO/01	F-104
36	Autotitrator	Metrohm	888 Titrande	CDTL/AUTO/03	GF-01
37	Autotitrator	Metrohm	888 Titrande	CDTL/AUTO/04	F-102
38	Autotitrator	Labindia	Titra smart	CDTL/AUTO/05	GF-01
39	Autotitrator	Labindia	Titra smart	CDTL/AUTO/06	GF-01
40	KF	Metrohm	795KFT	CDTL/KF/01	F-106
41	KF	Metrohm	890KFT	CDTL/KF/02	F-106
42	KF	Labndia	KAFI smart	CDTL/KF/03	GF-01
43	Coulometer	Metrohm	831KFC	CDTL/Coul/01	F-106
44	Atomic Absorption Spectrophotometer	GBC	Savant AA	CDTL/AAS/2	F-122
45	Autoclave	Osworld	NA	CDTL/MIC/AUTO/01	S-205
46	Autoclave	Medica Lab Equipment	325.91.6	CDTL/MIC/AUTO/02	S-225
47	Autoclave	Osworld	NA	CDTL/MIC/AUTO/03	S-225
48	Autoclave	Precision	PB-08	CDTL/AUTOCLAVE/04	S-225
49	Analytical Balance	Shimadzu	ATX 224	CDTL/BAL/03	F-106
50	Analytical Balance	LCGC Radwag	AS220/C/2	CDTL/BAL/04	F-104
51	Analytical Balance	LCGC Radwag	AS220/C/2	CDTL/BAL/06	F102
52	Analytical Balance	LCGC	377317/12	CDTL/BAL/07	S-202
53	Analytical Balance	Mettler toledo	AB204	CDTL/BAL/08	S-203
54	Analytical Balance	Mettler toledo	AT200	CDTL/BAL/09	S-204

55	Analytical Balance	Sartorius	BSA224S-CW	CDTL/BAL/11	S-224
56	Analytical Balance	Denver	SI-203	CDTL/BAL/10	S-220
57	Analytical Balance	Shimadzu	ATX 224	CDTL/BAL/12	G-04
58	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/15	F-128
59	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/16	F-106
60	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/17	S-204
61	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/18	G-02
62	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/19	S-204
63	Balance	Citizen	-----	CDTL/BAL/20	S-203
64	Balance	Citizen	-----	CDTL/BAL/21	F-102
65	Balance	Sartorius	Secura12510IN	CDTL/BAL/22	F-108
66	Balance	Sartorius	Secura12510IN	CDTL/BAL/23	F-108
67	BET Kinetic instrument	Lonza	ELX-808	CDTL/KIN/01	S-217
68	BOD Incubator	Osworld	3979	CDTL/BOD/01	S-223
69	BOD Incubator	Osworld	3980	CDTL/BOD/02	S-223
70	BOD Incubator	Metalab	JRIC-10	CDTL/BOD/03	S- 223
71	Biosafety Cabinet	-----	-----	CDTL/BSC/01	S-221
72	Break load test apparatus	Lab fine	-----	CDTL/COS-BLT/01	T-303
73	Centrifuge	Electrolab	CM-6MT	CDTL/CENTRI/01	F-101
74	Centrifuge	Remi motors	R-24	CDTL/CENTRI/02	G-010
75	Centrifuge	Remi motors	R-8C BL	CDTL/CENTRI/03	G-03
76	Centrifuge	Remi motors	R-8M PLUS	CDTL/CENTRI/04	F-105
77	Centrifuge	Remi motors	R-8C PLUS	CDTL/CENTRI/05	S-202
78	Centrifuge	Remi motors	R-8C PLUS	CDTL/CENTRI/06	T-304
79	Compressor	NA	NA	CDTL/Compsr/01	T-303
80	Dehumidifier	Bry-Air	NA	CDTL/DH/01	F-106
81	Dehumidifier	Bry-Air	NA	CDTL/DH/02	F-118
82	Disintegration apparatus	Veego	VTD-AV	CDTL/DT/02	S-202
83	Disintegration apparatus	Electrolab	EDI-2SA	CDTL/DT/05	F-119
84	Disintegration apparatus	Electrolab	EDI-2SA	CDTL/DT/06	T-304
85	Digital Force Gauge	Chattilon	DFA-010	CDTL/DFG/01	T-303
86	Digital Screw Gauge	Mitutoyo	293-561-30	CDTL/DSG/01	T-303
87	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/01	F-119
88	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/02	F-119
89	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/03	S-205
90	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/04	G-04
91	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/05	G-04
92	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/06	G-04
93	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/07	T-305
94	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/08	G-04
95	Dissolution apparatus	LabIndia	DS8001	CDTL/Disso/09	S-204
96	Dissoution media Prep sys	Lab inida	Media MET+	CDTL/DMP/01	G-04
97	Digital Flow meter for HPLC	GJC Instrument	5025000	CDTL/FM/01	F-125
98	Distillation unit(Food lab)	Velp Scientific	NA	NA	G-10
99	Drying time Machine	Labfine	NA	NA	T-304
100	Flame photometer	NA	-----	CDTL/FP/01	F-122
101	Foam Height Apparatus	NA	NA	CDTL/FHA/01	S-213
102	Freezer (Ultra deep)	SR lab Instruments	-----	CDTL/RF/039	S-217
103	FTIR	Shimadzu	IR Affinity - 1	CDTL/ /FTIR/02	F-118
104	FTIR	Perkin elmer	SP-BX-II	CDTL/ /FTIR/01	F-118
105	FTIR- Dry box	NA	-----	CDTL/FTIR-DB/01	F-118
106	FTIR-Press	NA	-----	CDTL/RP/01	118
107	FTIR-Press	NA	-----	CDTL/ /RP/02	118
108	FTIR-Press	Technosearch	-----	CDTL/ /RP/03	118
109	Fume Hood	Labguard	ST1500ABP	CDTL/FH/01	F-106
110	Fume Hood	Labguard	ST1500ABP	CDTL/FH/02	S-205
111	Fume Hood	Labguard	ST1500ABP	CDTL/FH/03	S-205
112	Fume Hood	Prama tech	NA	CDTL/FH/04	T-305
113	Funnel shaker for solvent extractor	Jiotech	ST1500ABP	CDTL/SFS/01	F-117

114	Gel Electrophoresis	Himedia	NA	CDTL/GEL/01	S-202
115	Heating mantle	NA	Na	CDTL/HM/01	S-204
116	Heating mantle	NA	NA	CDTL/HM/03	S-225
117	Heating mantle	NA	NA	CDTL/HM/04	F-107
118	Heating mantle	jindal	NA	CDTL/HM/05	S-205
119	Heating mantle	NA	NA	CDTL/HM/06	G-10
120	Heating mantle	Labquest	HME500	CDTL/HM/07	F-117
121	Heating mantle	Labquest	HME500	CDTL/HM/08	F-117
122	Heating mantle	Labquest	HME500	CDTL/HM/09	S-205
123	Heating mantle	Labquest	HME500	CDTL/HM/10	T-305
124	Heating Block	-----	NA	NA	S-217
125	Heating Block	-----	NA	NA	S-217
126	Hot plate digital	Tarson	NA	CDTL/HP/03	F-107
127	Hot plate digital	Tarson	NA	CDTL/HP/04	S-207
128	Incubator	Pathak		CDTL/INC/01	S-223
129	Incubator	Metlab	MSI-67	CDTL/INC-02	S-223
130	Incubator	Metlab	MSI-67	CDTL/INC-03	S-223
131	Laminar Air Flow	Air borne	NA	CDTL/LAF/01	S-219
132	Laminar Air Flow	Air borne	NA	CDTL/LAF/02	S-220
133	Laminar Air Flow	Air borne	NA	CDTL/LAF/03	S-221
134	Laminar Air Flow	Air borne	NA	CDTL/LAF/04	S-222
135	Magnetic stirrer	Emtec	NA	CDTL/MAG/02	S-202
136	Magnetic stirrer	Remi	NA	CDTL/MAG/03	G-04
137	Magnetic stirrer	Labquest	NA	CDTL/MAG/04	F-104
138	Magnetic stirrer	Labquest	NA	CDTL/MAG/05	T-305
139	Magnetic stirrer	Labquest	NA	CDTL/MAG/06	S-205
140	Melting point apparatus	LabIndia	MR-VIS	CDTL/MP/01	F-124
141	Microscope	Labomed	NA	CDTL/MIC/01	S-217
142	Microscope	Labomed	Vision2000	CDTL/MIC/02	S-217
143	Microscope	Labomed	Vision2001	CDTL/MIC/03	S-218
144	Mixer	Kay	NA	CDTL/MIX/02	T-305
145	Muffle furnance	NA	NA	CDTL/MF/01	205
146	Muffle furnance	NA	NA	CDTL/MF/02	205
147	Needle burner & Syringe destroyer	Hi Media	NA	CDTL/NSD/01	
148	Oil Bath	Metalab	NA	CDTL/OB/01	T-304
149	Oven	PATHAK	NA	CDTL/OV/02	F-107
150	Oven	METALAB	NA	CDTL/OV/01	F-101
151	Oven	JINDAL	NA	CDTL/OV/09	G-03
152	Oven	JINDAL	NA	CDTL/OV/10	G-04
153	Oven	MIC	NA	CDTL/OV/04	F-102
154	Oven	MIC	NA	CDTL/OV/08	S-202
155	Oven	Pathak	NA	CDTL/OV/05	S-205
156	Oven	Pathak	NA	CDTL/OV/06	S-224
157	Oven	MIC	NA	CDTL/OV/07	S-224
158	Oven	i therm	NA	CDTL/OV/03	S-205
159	Orbital shaker	Remi	NA	CDTL/SH/01	S-203
160	Orbital shaker	Remi	NA	CDTL/SH/02	G-10
161	Payoff measuring Tester	Labfine	NA	CDTL/PO/01	
162	pH Meter	LabIndia	Pico Plus	CDTL/pH/03	F-119
163	pH Meter	Euotech inst	pH 510	CDTL/pH/04	S-224
164	pH Meter	Euotech inst	Ph 700	CDTL/pH/09	F-104
165	pH Meter	ANALAB	pHcal10	CDTL/pH/12	F-104
166	pH Meter	Agilent	3200P	CDTL/pH/10	G-04
167	pH Meter	Agilent	3200P	CDTL/pH/11	S-202
168	pH Meter	Horiba	F72A	CDTL/pH/13	F-119
169	pH Meter	-----	-----	CDTL/pH/14	S- 201
170	Polarimeter	Rudolph research	-----	CDTL/POL/01	F-117
171	Polarimeter	Anton-Par	-----	CDTL/POL/02	F-117
172	Polarimeter	Anton-Par	-----	CDTL/POL/03	S-205

173	Profile projector	SIPCON	SP-300	CDTL/PP/02	T-304
174	Particle counter	-----	-----	CDTL/PC/01	S-203
175	Peel off Tester	Lab fine	-----	CDTL/COS-PT/01	T-303
176	Refractometer	Atago	RX-7000i	CDTL/RI/02	F-117
177	Scratch hardness tester	Labfine	NA	CDTL/SHT/01	ASN/Ramesh
178	Softening point apparatus	Labfine	NA	CDTL/SPA/01	ASN/Ramesh
179	Sonicattor	Pci service	NA	CDTL/Soni/02	F-104
180	Sonicattor	MIC	-----	CDTL//Soni/04	T-304
181	Sonicattor	MIC	-----	CDTL//Soni/05	S-203
182	Sonicattor	Labdevice	NA	CDTL//Soni/07	F101
183	Sonicattor	MIC	-----	CDTL//Soni/08	S-202
184	Sonicattor	MIC	-----	CDTL//Soni/09	S-204
185	Sonicattor	Pci service	NA	CDTL//Soni/10	G-03
186	Sonicattor	Pci service	NA	CDTL//Soni/11	F-101
187	Solvent cabinet	Prama	-----	-----	G-011
188	Solvent cabinet	Prama	-----	-----	G-011
189	Spot Extractor	Prama	-----	CDTL/SPEX/01	F-104
190	Spot Extractor	Prama	-----	CDTL/SPEX/02	F-104
191	Spot Extractor	Prama	-----	CDTL/SPEX/03	F-101
192	Spot Extractor	Prama	-----	CDTL/SPEX/04	F-101
193	Spot Extractor	Prama	-----	CDTL/SPEX/05	S-202
194	Spot Extractor	Prama	-----	CDTL/SPEX/06	S-203
195	Spot Extractor	Prama	-----	CDTL/SPEX/07	S-204
196	Standard Weight Box	IDEMI Mumbai	NA	-----	QA
197	Tensile compression Device	Chattilon	TCD200	CDTL//TCD/01	T-303
198	Tensile compression Device	Chattilon	TCD200	CDTL//TCD/02	T-303
199	Tensile compression Device	Mecmesin	AFG-100N	CDTL//TCD/03	T-303
200	Traceble Digital Calliper	Fisher Scientific	NA	CDTL//TDC/01	T-303
201	Traceble Digital Timer	Fisher Scientific	NA	CDTL//TDT/01	T-303
202	Tyvec pouch Burst strength	NA	NA	CDTL/TPBST/01	S- 203
203	TOC Analyzer	shimazu	TOC-LCPH	CDTL//TOC/1	G-02
204	Tablet Friability Tester	lab india	-----	CDTL//TFT/01	G-013
205	Tablet Hardnes Tester	lab india	-----	CDTL//THT/01	G-013
206	UV/VIS SPECTROPHOTOMETER	LabIndia	UV300+	CDTL//UV/01	G-03
207	UV/VIS SPECTROPHOTOMETER	LabIndia	UV300+	CDTL//UV/02	S-203
208	UV/VIS SPECTROPHOTOMETER	Perkin Elmer	UV150-02	CDTL//UV/03	F-117
209	Ultra violet Nail paint Dryer	-----	-----	CDTL/COS-ND/01	T-303
210	Vaccume oven	Osworld	NA	CDTL//VO/01	F-105
211	Vaccume oven	LAB HOSP	NA	CDTL//VO/02	F-107
212	Vaccume oven	HMG India	NA	CDTL//VO/03	F-105
213	Vaccume oven	HMG	NA	CDTL//VO/04	S-203
214	Vaccume oven	Fourtech	NA	CDTL//VO/05	S-206
215	Vaccume oven	SM Scientific	NA	CDTL//VO/06	F-105
216	Vaccume oven	META LAB	NA	CDTL//VO/07	G-02
217	Vaccume oven	JEIOTECH	ov-12	CDTL//VO/08	F-105
218	Vortex	SPINIX	NA	CDTL//VOR/01	F-102
219	Vortex	SPINIX	NA	CDTL//VOR/02	F-202
220	Vortex	Remi	NA	CDTL//VOR/03	S-204
221	Vortex	SPINIX	NA	CDTL//VOR/04	S-220
222	Vortex	Remi	NA	CDTL//VOR/05	S-220
223	Vortex	-----	-----	CDTL//VOR/07	S-217
224	Vortex	Remi	CM101	CDTL//VOR/08	F104
225	Vortex	Remi	CM102	CDTL//VOR/09	G-03
226	Vortex	Labquest	VMSP	CDTL//VOR/10	F104
227	Vortex	Labquest	VMSP	CDTL//VOR/11	T-305
228	Vortex	Labquest	VMSP	CDTL//VOR/12	S-205
229	Vortex	Labquest	VMSP	CDTL//VOR/13	G-10
230	Vortex mix	Remi	NA	CDTL//VOR/14	G-04
231	Water Bath-Shaker	-----	NA	CDTL//SWB/01	S-202

232	Water Bath-Shaker	-----	NA	CDTL/SWB/02	F-107
233	Water Bath	-----	NA	CDTL/WB/06	S-224
234	Water Bath	-----	NA	CDTL/WB/10	S-205
235	Water Bath	-----	NA	CDTL/WB/11	F-107
236	Water Bath	Foodlab	NA	CDTL/WB/12	F-107
237	Water Bath	Foodlab	NA	CDTL/WB/13	T-305
238	Water Bath	Labquest	WBC-012	CDTL/WB/14	F-107
239	Water Bath	Labquest	WBC-012	CDTL/WB/15	S-205
240	Water Bath	Labquest	WBC-012	CDTL/WB/16	T-305
241	Water Bath	Labquest	WBC-012	CDTL/WB/17	G-10
242	Water purifier	Pure Lab Ultra	ELGA DV-35	CDTL//WP/01	T-304
243	Water purifier	Millipore	Q815 V	CDTL/WP/02	F-101
244	Water purifier	Millipore	Elix	CDTL/WP/05	F-105
245	Water purifier	Thermo	Genpure	CDTL/WP/04	G-13
246	Water purifier	Thermo	Genpure	CDTL/WP/03	S-211
247	Zone Reader	-----	-----	CDTL//ZR/01	S-227